

HOTLINE: Effective May 20, 2019

New Test	3001479	Respiratory Viral Panel by PCR	RVPPCR
-----------------	----------------	---------------------------------------	---------------

Methodology: Qualitative Polymerase Chain Reaction

Performed: Sun-Sat

Reported: 1-3 days

Specimen Required: Collect: Nasopharyngeal Swab.

Specimen Preparation: Place in viral transport media (ARUP Supply #12884). Available online through eSupply using ARUP Connect[™] or contact ARUP Client Services at (800) 522-2787. Place each specimen in an individually sealed bag.

Storage/Transport Temperature: Frozen

Unacceptable Conditions: Specimens not in viral transport media.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 4 days; Frozen: 1 month

Interpretive Data: A negative result does not rule out the presence of PCR inhibitors in the patient specimen or assay-specific nucleic acid in concentrations below the level of detection by this assay.

Note: This test detects influenza A, influenza B, RSV, human metapneumovirus, human rhinovirus, and adenovirus. Detects and differentiates parainfluenza 1, 2, 3, and 4.

CPT Code(s): 87632

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.