

New Test

3001535

TPMT and NUDT15

TPMT2



Additional Technical Information



Supplemental Resources

Methodology: Polymerase Chain Reaction/Fluorescence Monitoring
Performed: Varies
Reported: 5-10 days

Specimen Required: Collect: **Whole Blood:** Lavender (EDTA), Pink (K₂EDTA), or Yellow (ACD Solution A or B).
Saliva: Collection Device by DNA Genotek (OCD-100, ARUP Supply #49295) available online through eSupply using ARUP Connect™ or by contacting ARUP Client Services at (800) 522-2787.
Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL) OR Transport the Saliva Collection Device.
Storage/Transport Temperature: **Whole Blood:** Refrigerated.
Saliva: Room temperature.
Unacceptable Conditions: Plasma or serum. Specimens collected in sodium heparin or lithium heparin.
Stability (collection to initiation of testing): **Whole Blood:** Ambient: 72 hours; Refrigerated: 1 week; Frozen: 1 month
Saliva: Ambient: 2 weeks; Refrigerated: Unacceptable; Frozen: Unacceptable

Reference Interval: By report

Interpretive Data:

Background Information for TPMT and NUDT15:

Characteristics: Thiopurine drug therapy is used for autoimmune diseases, inflammatory bowel disease, acute lymphoblastic leukemia, and to prevent rejection after solid organ transplant. The inactivation of thiopurine drugs is catalyzed in part by thiopurine methyltransferase (TPMT) and nudix hydrolase 15 (NUDT15). Variants in the *TPMT* and/or *NUDT15* genes are associated with an accumulation of cytotoxic metabolites leading to increased risk of drug-related toxicity with standard doses of thiopurine drugs. These effects on thiopurine catabolism can be additive.

Inheritance: Autosomal co-dominant.

Cause: *TPMT* and *NUDT15* variants affect enzyme expression or activity.

Variants Tested: See the "Additional Technical Information" document.

Clinical Sensitivity: 95 percent.

Methodology: Polymerase chain reaction (PCR) and fluorescence monitoring.

Analytical Sensitivity and Specificity: 99 percent.

Limitations: Only the targeted *TPMT* and *NUDT15* variants will be detected by this test. Because the complex TPMT*3A allele contains the variants found in the *3B and *3C alleles, this test cannot distinguish the 3A/Negative genotype (intermediate enzyme activity) from the rare *3B/*3C genotype (no or low enzyme activity). Genotyping may reflect donor status in patients who have received allogenic stem cell or bone marrow transplants within 2 weeks of specimen collection. Actual enzyme activity and expression and risk for adverse reactions to thiopurines may be affected by additional genetic and non-genetic factors not evaluated by this test. Diagnostic errors can occur due to rare sequence variations. Genotyping does not replace the need for therapeutic drug monitoring and clinical observation.

See Compliance Statement C: www.aruplab.com/CS

Note: Whole blood is the preferred specimen. Saliva samples that yield inadequate DNA quality and/or quantity will be reported as inconclusive if test performance does not meet laboratory-determined criteria for reporting.

CPT Code(s): 81335; 81306

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.