

HOTLINE: Effective August 19, 2019

New Test 3001760 Pancreatitis (PRSS1) Deletion/Duplication

PRSS1 DD



Patient History for Pancreatitis (*PRSS1*) Del/Dup



Additional Technical Information

Methodology: Multiplex Ligation-dependent Probe Amplification

Performed: Varies **Reported:** 12-14 days

Specimen Required: Collect: Lavender (K₂EDTA), Pink (K₂EDTA), or Yellow (ACD Solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 2 mL)

Storage/Transport Temperature: Refrigerated.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 6 months

Reference Interval: By Report

Interpretive Data:

Background Information for Pancreatitis (PRSS1) Deletion/Duplication:

Characteristics: Hereditary pancreatitis typically presents in late childhood with recurrent episodes of pancreatic inflammation, abdominal pain, nausea, and vomiting. Ultimately, this evolves into chronic pancreatitis resulting in permanent pancreatic damage.

Epidemiology: Incidence of chronic pancreatitis is 5-12 in 100,000 per year and prevalence is approximately 50 in 100,000.

Inheritance of SPINK1-related pancreatitis: Autosomal recessive and possibly digenic.

Penetrance: Variable.

Cause: Pathogenic variants in PRSS1, SPINK1, CFTR, CASR, CTRC, CPA1 and CLDN2 genes are associated with pancreatitis.

Clinical Sensitivity: Unknown.

Methodology: Multiplex ligation-dependent probe amplification (MLPA) of the PRSS1 gene.

Analytical Sensitivity and Specificity: 99 percent.

Limitations: Diagnostic errors can occur due to rare sequence variations. Single base pair substitutions, small deletions/duplications, regulatory region and deep intronic variants will not be detected. Deletion/duplication breakpoints will not be determined. Variants in genes other than *PRSS1* will not be detected.

Counseling and informed consent are recommended for genetic testing. Consent forms are available online at www.aruplab.com See Compliance Statement C: www.aruplab.com/CS

New York DOH approval pending. Call for status update.

CPT Code(s): 81479

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.