

New Test **3001783** **Dermatomyositis and Polymyositis Panel** **COMBI PAN**



Additional Technical Information

Methodology: Qualitative Immunoprecipitation/Semi-Quantitative Multiplex Bead Assay/Qualitative Immunoblot
Performed: Sun-Sat
Reported: 7-18 days

Specimen Required: Collect: Serum Separator Tube (SST).
Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer two 1 mL serum aliquots to ARUP Standard Transport Tubes. (Min: 0.5 mL/aliquot)
Storage/Transport Temperature: Refrigerated.
Unacceptable Conditions: Hemolyzed, hyperlipemic, icteric, heat-treated or contaminated specimens.
Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

Reference Interval:

Test Number	Components	Reference Interval	
0099592	Jo-1 Antibody, IgG	29 AU/mL or less	Negative
		30-40 AU/mL	Equivocal
		41 AU/mL or greater	Positive
	PL-7 (threonyl-tRNA synthetase) Antibody	Negative	
	PL-12 (alanyl-tRNA synthetase) Antibody	Negative	
	EJ (glycyl-tRNA synthetase) Antibody	Negative	
	SRP (Signal Recognition Particle) Ab	Negative	
	OJ (isoleucyl-tRNA synthetase) Antibody	Negative	
	Mi-2 (nuclear helicase protein) Antibody	Negative	
	P155/140 Antibody	Negative	
	SAE1 (SUMO activating enzyme) Ab	Negative	
	MDA5 (CADM-140) Ab	Negative	
	NXP2 (Nuclear matrix protein-2) Ab	Negative	
	TIF1-gamma Ab	Negative	

Interpretive Data: Refer to report.
 See Compliance Statement D: www.aruplab.com/CS

Note: Antibodies: PL-7, PL12, EJ, OJ, SRP, Jo-1, Mi-2, P155/140, SAE1, MDA5, NXP2, TIF1-gamma

CPT Code(s): 83516 x11; 86235

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.