

HOTLINE: Effective August 19, 2019

New Test **3001831** **Troponin T (cTnT) 5th Generation** **CTNT**

Methodology: Quantitative Electrochemiluminescent Immunoassay (ECLIA)
Performed: Mon, Wed, Fri
Reported: 1-4 days

Specimen Required: Collect: Plasma Separator Tube (PST), or Green (Lithium Heparin).
Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: **CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.**
Unacceptable Conditions: Specimens collected in potassium oxalate, sodium fluoride, or sodium citrate. Grossly hemolyzed specimens.
Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 24 hours; Frozen: 1 year

Reference Interval:

| | |
|--------|-------------------------------|
| Female | Less than or equal to 10 ng/L |
| Male | Less than or equal to 15 ng/L |

CPT Code(s): 84484

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.