

TEST CHANGE

Chronic Lymphocytic Leukemia Mutation Panel by Next Generation Sequencing

3001858, CLL NGS	
Specimen Requirements:	
Patient Preparation:	
Collect:	Lavender (EDTA), Green (sodium heparin), Bone Marrow (EDTA), or Bone Marrow (sodium heparin). Fresh-frozen tissue. <u>New York State Clients: Lavender (EDTA)</u>
Specimen Preparation:	 Whole Blood and Bone Marrow: Transport 3 mL. (Min: 1.5 mL) Fresh-frozen Tissue: Transport 5 mg fresh-frozen tissue. (Min: 5 mg) Separate specimens must be submitted when multiple tests are ordered. <u>New York State Clients: Transport 5 mL whole blood (Min: 2 mL) or 2 mL bone marrow (Min: 2 mL).</u>
Transport Temperature:	Whole Blood or Bone Marrow: Refrigerated. Fresh-frozen Tissue: Frozen.
Unacceptable Conditions:	Serum, plasma, grossly hemolyzed specimens, buccal brush or swab, FFPE tissue.
Remarks:	
Stability:	Whole Blood or Bone Marrow: Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable Fresh-frozen Tissue: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month <u>New York State Clients: Ambient: 72 hours; Refrigerated: 4</u> days; Frozen: Unacceptable
Methodology:	Massively Parallel Sequencing
Performed:	Varies
Reported:	12-14 days
Note:	Genes tested: ATM; BCL2; BIRC3*; BRAF; BTG1; BTK; CARD11; CD79B; CXCR4; DDX3X; FBXW7; IKZF3; KRAS; MAP2K1; MED12; MGA; MYD88; NOTCH1; NRAS; PLCG2; POT1; RNASEH2A; RNASEH2B; RPS15*; SAMHD1; SF3B1; TP53; XPO1; ZMYM3 *One or more exons are not covered by sequencing for the indicated gene; see Additional Technical Information test fact sheet.
CPT Codes:	81450



New York DOH Approval Status:

Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Refer to report.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA-certified laboratory and is intended for clinical purposes.

Reference Interval:

By report