

HOTLINE: Effective November 18, 2019

---

**New Test**      **3002008**      **Human Papillomavirus (HPV) High Risk by in situ Hybridization, Paraffin**      **HPVHR ISH**

---

**Methodology:** In situ Hybridization

**Performed:** Mon-Fri

**Reported:** 2-5 days

**Specimen Required:** Collect: Tissue.

Specimen Preparation: Formalin fix (10 percent neutral buffered formalin) and paraffin-embed tissue. Transport tissue block or 5 unstained 5-micron slides in a tissue transport kit (recommended but not required) (ARUP supply #47808). Available online through eSupply using ARUP Connect<sup>™</sup> or contact ARUP Client Services at (800) 522-2787. (Min: 4 slides) Protect paraffin block and/or slides from excessive heat.

Storage/Transport Temperature: Room temperature or refrigerated. Ship in cooled container during summer months.

Remarks: Include surgical pathology report.

Unacceptable Conditions: Specimens fixed or processed in alternative fixatives (alcohol, Prefer) or heavy metal fixatives (B-4 or B-5). Frozen specimens.

Stability (collection to initiation of testing): Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

**Interpretive Data:**

Refer to report.

See Compliance Statement B: [www.aruplab.com/CS](http://www.aruplab.com/CS)

**CPT Code(s):** 88365

New York DOH Approved.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.