

TEST CHANGE

Anti-gp210 Antibody, IgG 3002477, GP210 AB

Specimen Requirements:

Patient Preparation:

Collect: Serum <u>separator tube</u> (SST).

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection.

Transfer 0.5 mL serum to ARUP <u>standard transport</u>

tube. Standard Transport Tube. (Min: 0.1 mL)

Transport Temperature: Refrigerated

Unacceptable Conditions: Nonserum Non-serum, heat-inactivated, contaminated, grossly

icteric, severely lipemic, grossly hemolyzed specimens, or

Effective Date: July 21, 2025

inclusion of fibrin clot.

Remarks:

Stability: After separation from cells: Ambient: 48 hours; Refrigerated: 2

weeks; Frozen: 1 month (avoid repeated freeze/thaw cycles)

Methodology: Semi-Quantitative Enzyme-Linked Immunosorbent Assay

Performed: <u>Thu</u>Wed

Reported: 1-8 days

Note:

CPT Codes: 83516

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

GP210 IgG antibodies can be detected in patients with primary biliary cholangitis (PBC) and may be of diagnostic relevance in a subset of patients with PBC who are negative for anti-mitochondrial antibodies (AMA). These antibodies have a relatively low sensitivity with excellent specificity for PBC. A negative result does not rule out PBC.

Component Interpretation

Anti-gp210 20.1 Units or less

Antibody, IgG Negative 20.124.9 Units

24.9 Units Equivocal 25.0 Units or greater Positive

Reference Interval:



Effective Date: July 21, 2025

Test Number	Components	Reference Interval
	Anti-gp210 Antibody, IgG	24.9 Units or less