

HOTLINE: Effective May 18, 2020

New Test	3002568	Monoclonal Protein Study, Serum	IFE SPEP
Methodology: Performed: Reported:	Qualitative Immu Sun-Sat 1-5 days	nofixation Electrophoresis/Quantitative Capillary Electrophoresis/Quantitative Spectrophotomet	ry

Specimen Required: Collect: Serum Separator Tube (SST).

 Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1.5 mL serum to an ARUP Standard Transport Tube. (Min: 1 mL)

 Storage/Transport Temperature: Refrigerated.

 Unacceptable Conditions: Plasma.

 Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 1 month

Reference Interval:

Test Number	Components	Reference Interval		
0050640	Protein Electrophoresis, Serum	Effective August 19, 2019		
		Test Number	Components	Reference Interval
			Total Protein, Serum	Refer to report
			Albumin	Refer to report
			Alpha-1 Globulins	Refer to report
			Alpha-2 Globulins	Refer to report
			Beta Globulins	Refer to report
			Gamma	Refer to report

Note: A copy of the graph will follow the final report.

CPT Code(s): 84155; 84165; 86334

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.