

HOTLINE: Effective May 17, 2021

**New Test**      **3002581**      **Vaginitis Panel by TMA**      **VPAN TMA**

**Methodology:** Qualitative Transcription-Mediated Amplification  
**Performed:** Tue, Thu, Sat  
**Reported:** 1-4 days

**Specimen Required:** Patient Prep: Patient must be 14 years of age or older.  
Collect: Vaginal specimen collected with pink swab from Aptima MultiTest Swab Collection kit (ARUP supply #55224 PK/50 or #55229 PK/10) available online through eSupply using ARUP Connect or contact Client Services at (800) 522-2787.  
Specimen Preparation: Place swab in MultiTest Swab Specimen Transport Tube, break shaft at scoreline, then recap tube.  
Storage/Transport Temperature: Refrigerated.  
Unacceptable Conditions: Specimens in any transport media other than indicated above. Specimen in MultiTest swab transport media without a swab.  
Stability (collection to initiation of testing): Ambient: 30 days; Refrigerated: 30 days; Frozen: 90 days

**Reference Interval:**

Test Number	Components	Reference Interval
3002582	Bacterial Vaginosis by Transcription-mediated Amplification (TMA)	Negative
3002583	Candida glabrata, Candida species, and Trichomonas vaginalis by Transcription-mediated Amplification (TMA)	Negative

**Interpretive Data:**

See report

**CPT Code(s):** 81513, 87481, 87661

New York DOH Approved.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.