

**NEW TEST** 

**Click for Pricing** 

Hereditary Melanoma Panel, Sequencing and Deletion/Duplication

3002673, MELCAN NGS

Specimen Requirements:

Patient Preparation:

Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 2 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens;

saliva; buccal brush, or swab; FFPE tissue.

Remarks:

Stability: Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

Methodology: Massively Parallel Sequencing/Sequencing

Performed: Varies

Reported: 3 weeks

Note: GENES TESTED: BAP1; BRCA2; CDK4; CDKN2A\*; MC1R; MITF\*;

POT1; PTEN\*; RB1\*; TERT; TP53 \*One or more exons are not covered by sequencing and/or deletion/duplication analysis for the indicated gene; see Additional Technical Information.

Effective Date: February 21, 2023

CPT Codes: 81167; 81216; 81404; 81321; 81323; 81351; 81479

New York DOH Approval Status: Specimens from New York clients will be sent out to a New

York DOH approved laboratory, if possible.

Interpretive Data:

Refer to report.

Reference Interval:

By report

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.