

HOTLINE: Effective November 16, 2020

New Test **3002976** **Adenosine Deaminase, Serum or Plasma** **ADASP**

Methodology: Quantitative Spectrophotometry
Performed: Sun, Tue, Thu
Reported: 1-4 days

Specimen Required: Collect: Plain red, SST, or Green Heparin Plasma (Lithium or Sodium)
Specimen Preparation: Allow serum tube to clot completely at room temperature. Separate serum or plasma from cells within 4 hours of collection. Transfer 0.5 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.2 mL)
Storage/Transport Temperature: Frozen. Specimen must remain frozen until received in lab.
Remarks: Indicate source on requisition. Conditions: EDTA, citrate or oxalate.
Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 1 week; Frozen: 1 month

Reference Interval: 0-15 U/L

Interpretive Data:

See Compliance Statement B : www.aruplab.com/CS

CPT Code(s): 84311

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.