

HOTLINE: Effective November 16, 2020

New Test **3002980** **Adenosine Deaminase in Peritoneal Fluid** **ADAPERIT**

Methodology: Quantitative Spectrophotometry
Performed: Sun, Tue, Thu
Reported: 1-4 days

Specimen Required: Patient Prep: Collect specimens in leak proof container.
Collect: Peritoneal Fluid.
Specimen Preparation: Centrifuge specimen at room temperature. Transfer 0.5mL fluid to an ARUP Standard Transport Tube and freeze. (Min: 0.2 mL)
Storage/Transport Temperature: Frozen. Specimen must remain frozen until received in lab.
Remarks: Indicate source on requisition.
Unacceptable Conditions: Whole blood. Bronchoalveolar Lavage (BAL) specimens. Turbid specimens.
Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 1 week; Frozen: 1 month

Reference Interval: 0-30 U/L

Interpretive Data:

See Compliance Statement B www.aruplab.com/CS

CPT Code(s): 84311

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.