

HOTLINE: Effective November 16, 2020

New Test 3003017 Autoimmune Neuromuscular Junction Reflexive Panel MUWA R2



Additional Technical Information



Supplemental Resources

Methodology: Quantitative Radioimmunoassay/Qualitative Radiobinding Assay/Semi-Quantitative Flow Cytometry/Semi-Quantitative Indirect

Fluorescent Antibody

Performed: Refer to individual components

Reported: 2-8 days

Specimen Required: Collect: Serum separator tube (SST).

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transport 2 mL serum. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

<u>Unacceptable Conditions:</u> Plasma. Contaminated, hemolyzed, or severely lipemic specimens. <u>Stability (collection to initiation of testing):</u> Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

Reference Interval:

| Test Number | Components | Reference Interva | al |
|-------------|--|-------------------|----------------------------------|
| 0080009 | Acetylcholine Receptor Binding Antibody | | |
| | | Negative | 0.0-0.4 nmol/L |
| | | Positive | 0.5 nmol/L or greater |
| | | | |
| 0099580 | Acetylcholine Receptor Blocking Antibody | | |
| | | Negative: | 0-26 percent blocking |
| | | Indeterminate: | 27-41 percent blocking |
| | | Positive: | 42 percent or greater blocking |
| 0000521 | A control de lieu De contro Medicine Antibodo | | |
| 0099521 | Acetylcholine Receptor Modulating Antibody | Negative | 0-45 percent modulating |
| | | | |
| | | Positive | 46 percent or greater modulating |
| | | | |
| 0092628 | P/Q-Type Voltage-Gated Calcium Channel (VGCC) Antibody | | |
| | | Negative | 0.0 to 24.5 pmol/L |
| | | Indeterminate | 24.6 to 45.6 pmol/L |
| | | Positive | 45.7 pmol/L or greater |
| | N-Type Voltage-Gated Calcium Channel (VGKC) Antibody | | |
| | | Negative | 0.0 to 69.9 pmol/L |
| | | Indeterminate | 70.0 to 110.0 pmol/L |
| | | Positive | 110.1 pmol/L or greater |
| 3003020 | Ganglionic Acetylcholine Receptor Antibody | | |
| | | Negative | pmol/L or less |
| | | Indeterminate | pmol/L |
| | | Positive | pmol/L or greater |
| 2004890 | Voltage-Gated Potassium Channel (VGKC) Antibody, Serum | | |
| 2004890 | vonage-Galeu Fotassium Chaimei (VGKC) Antibody, Sefum | Negative | 31 pmol/L or less |
| | | Indeterminate | 32-87 pmol/L |
| | | Positive | 88 pmol/L or greater |
| | | | <u> </u> |
| 2005636 | Titin Antibody | | |
| | | Negative | 0.00-0.45 IV |
| | | Indeterminate | 0.46-0.71 IV |
| | | Positive | 0.72 IV or greater |
| 0050746 | Striated Muscle Antibodies, IgG with Reflex to Titer | Less than 1:40 | |



HOTLINE: Effective November 16, 2020

| 2009456 | Leucine-Rich, Glioma-Inactivated Protein 1 Antibody, IgG with Reflex to Titer, Serum | Less than 1:10 |
|---------|--|----------------|
| 2009452 | Contactin-Associated Protein-2 Antibody, IgG with Reflex to Titer, Serum | Less than 1:10 |

Interpretive Data:

Refer to report.

See Compliance Statement A: www.aruplab.com/CS

Note: If Acetylcholine Receptor Binding Antibody result is greater than 0.4 nmol/L or Acetylcholine Receptor Blocking Antibody result is greater than 26 percent, then Acetylcholine Receptor Modulating Antibody will be added. If Striated Muscle Ab is detected, then a titer will be added. If VGKC is Indeterminate or Positive, LGI1 Antibody IgG and CASPR2 Antibody IgG will be added. If LGI1 antibody IgG is positive, then LGI1 antibody IgG titer will be added. If CASPR2 antibody IgG is positive, then LGI1 antibody IgG titer will be added. Additional charges apply.

CPT Code(s): 83519 x5; 83516 x2; 86255; if reflexed, add 83516 and/or 86256 and/or 86255x2, if further reflexed add 86256 per titer

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.