

New Test **3003477** **Membranous Nephropathy Comprehensive Autoantibody Panel** **MNCA PAN**
 Available Now



Additional Technical Information

Methodology: Semi-Quantitative Indirect Fluorescent Antibody
Performed: Tue
Reported: 1-8 days

Specimen Required: Collect: Serum Separator Tube
Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.2 mL)
Storage/Transport Temperature: Refrigerated
Unacceptable Conditions: Hemolyzed, hyperlipemic, icteric, heat-treated or contaminated
Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month

Reference Interval:

2011828	Phospholipase A2 Receptor (PLA2R) Antibody, IgG with Reflex to Titer	Less than 1:10
3003480	Antithrombospondin Type-1 Domain-Containing 7A (THSD7A) Antibody, IgG with Reflex to Titer	Less than 1:10

Interpretive Data:

Refer to individual components.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Note: If Phospholipase A2 Receptor Antibody, IgG is positive, then a Phospholipase Receptor A2 Antibody, IgG titer will be added. Additional charges apply. If THSD7A Antibody, IgG is positive, then a THSD7A Antibody, IgG titer will be added.. Additional charges apply.

CPT Code(s): 86255 x2; if reflexed, add 86256 x2

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.