

HOTLINE: Effective February 16, 2021

---

<b>New Test</b>	<b>3003525</b>	<b>Kappa/Lambda Light Chain Panel by In Situ Hybridization Stain Only</b>	<b>SO KLP ISH</b>
-----------------	----------------	---	-------------------

**Methodology:** In situ hybridization (ISH)  
**Performed:** Mon-Fri  
**Reported:** 2-5 days

**Specimen Required:** Collect: Tissue or cells.

Specimen Preparation: Formalin fix (10 percent neutral buffered formalin) and paraffin-embed specimen (cells must be prepared into a cellblock). Protect paraffin block and/or slides from excessive heat. Transport tissue block or 8 unstained (3- to 5-micron thick sections), positively-charged slides in a tissue transport kit (ARUP supply #47808). Available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787. (Min: 4 slides). If sending precut slides, do not oven bake.

Storage/Transport Temperature: Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

Remarks: **IMMUNOHISTOCHEMISTRY ORDERING AND SUBMISSION DETAILS:** Submit electronic request. If you do not have electronic ordering capability, use an ARUP Immunohistochemistry Stain Form (#32978) with an ARUP client number. For additional technical details, contact ARUP Client Services at (800) 522-2787.

Unacceptable Conditions: **Unacceptable Conditions:** Specimens submitted with non-representative tissue type. Depleted specimens.

Stability (collection to initiation of testing): Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

**Interpretive Data:**

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

**CPT Code(s):** 88368; 88369

New York DOH Approved.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.