

HOTLINE: Effective June 1, 2021

New Test **3003648** **COVID-19 IgG, Semi-Quantitative by CIA** **COV19G SQ**
 Available Now

Methodology: Semi-Quantitative Chemiluminescent Immunoassay
Performed: Sun-Sat
Reported: Within 24 hours

Specimen Required: Collect: Serum separator tube (SST). Also acceptable: EDTA or lithium heparin.
Specimen Preparation: Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: Refrigerated.
Unacceptable Conditions: Specimens containing particulate material or otherwise obviously contaminated. Severely hemolyzed, heat-inactivated, or lipemic specimens.
Stability (collection to initiation of testing): After separation from cells: Ambient: 72 hours; Refrigerated: 1 week; Frozen: 1 month (avoid repeated freeze/thaw cycles)

Reference Interval:

Less than 1.00 Index Value	Negative
Greater than or equal to 1.00 Index Value	Positive

Interpretive Data:

This test is for in vitro diagnostic use under the FDA Emergency Use Authorization (EUA). In compliance with this authorization, please visit <https://www.aruplab.com/infectious-disease/coronavirus/testing> for more information and to access the applicable fact sheets.

CPT Code(s): 86769

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.