

TEST CHANGE

Propoxyphene and Metabolite, Urine

3003726, PROPOX U

Specimen Requirements:

Patient Preparation:

Collect: Urine

Specimen Preparation: Transfer 2 mL urine to an ARUP **standard transport tube**~~Standard Transport Tube~~. (Min: 0.7 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: **Frozen**~~Refrigerated~~. Also acceptable: **Refrigerated**~~Frozen~~

Unacceptable Conditions:

Remarks:

Stability: Ambient: **3 days**~~72 hours~~; Refrigerated: 1 week; Frozen: 3 months

Methodology: Quantitative Gas Chromatography-Mass Spectrometry (GC-MS)

Note: Amitriptyline is a known interference.

CPT Codes: 80367 ~~(Alt code: G0480)~~

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer to~~By~~ report