

HOTLINE: Effective August 16, 2021

New Test	3003801JAK2 (V617F) Mutation by ddPCR, Qualitative with Reflex to JAK2 Exon 12 Mutation Analysis by PCR	PV RFX
ii,	Additional Technical Information	
Methodology:		
Performed:	DNA Isolation: Sun-Sat Assay: Varies	
Reported:	3-12 days	
Specimen Requir	ired: Collect: Whole blood or bone marrow: Lavender (EDTA), preferred. Also acceptable: Green (sodium heparin).	
	Specimen Preparation: Whole Blood: Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)	
	Bone Marrow: Do not freeze. Transport 3 mL bone marrow. (Min: 1 mL)	
	Storage/Transport Temperature: Refrigerated.	
	<u>Unacceptable Conditions:</u> Plasma, serum, FFPE tissue blocks/slides, or fresh or frozen tissue, DNA extracted by a no	
	lab. Specimens collected in anticoagulants other than EDTA or sodium heparin. Clotted or grossly hemolyzed specim Stability (collection to initiation of testing): Pefrigameted: 7 days: Frazer, Unsecontable	nens.
	Stability (collection to initiation of testing): Refrigerated: 7 days; Frozen: Unacceptable	

Interpretive Data:

Refer to report.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Note: If JAK2 QUAL is reported as "Not Detected" then JAK2 Exon 12 Mutation Analysis will be added. Additional charges apply.

CPT Code(s): 81270; if reflexed, add 81279

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.