

HOTLINE: Effective February 22, 2022

---

**New Test**      **3004383**      **Cerebral Autosomal Dominant Arteriopathy with Subcortical Infarcts and Leukoencephalopathy, CADASIL (*NOTCH3*), Sequencing**      **NOTCH3 NGS**



Additional Technical Information



Patient History for CADASIL (*NOTCH3* Gene) Testing

**Methodology:** Massively Parallel Sequencing  
**Performed:** Varies  
**Reported:** 3 weeks

**Specimen Required:** Collect: Lavender (EDTA) or Yellow (ACD Solution A or B).  
Specimen Preparation: Transport 5 mL whole blood. (Min: 3 mL)  
Storage/Transport Temperature: Refrigerated.  
Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens  
Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

**Reference Interval:** By report

**Interpretive Data:**  
Refer to report.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA-certified laboratory and is intended for clinical purposes.

Counseling and informed consent are recommended for genetic testing. Consent forms are available online.

**Note:** Gene Tested: *NOTCH3* (NM\_000435)  
Exon 1 is not covered by sequencing.

**CPT Code(s):** 81406

New York DOH approval pending. Call for status update.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.