

**TEST CHANGE**

**Paraneoplastic Reflexive Panel, CSF**

3004517, PNSPAN CSF

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** CSF.

**Specimen Preparation:** Transfer 2 mL CSF to an ARUP standard transport tube. (Min: 1 mL).

**Transport Temperature:** Refrigerated

**Unacceptable Conditions:** Contaminated, heat-inactivated, hemolyzed, or lipemic specimens

**Remarks:**

**Stability:** Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month

**Methodology:** Semi-Quantitative Indirect Fluorescent Antibody (IFA)/Qualitative Immunoblot/Semi-Quantitative Cell-Based Indirect Fluorescent Antibody

**Performed:** Wed

**Reported:** 1-9 days

**Note:** Purkinje Cell (PCCA) antibody and ~~neuronal nuclear~~ **Neuronal Nuclear** (ANNA) antibody IgG are screened by IFA. If the IFA screen is indeterminate, then a Neuronal Nuclear Antibodies (Hu, Ri, Yo, Tr/DNER) IgG by Immunoblot will be added. If the IFA screen is positive at 1:1 or greater, then a PCCA/ANNA antibodies titer and Neuronal Nuclear Antibodies (Hu, Ri, Yo, Tr/DNER) IgG by Immunoblot will be added. Additional charges apply. If CV2 Antibody IgG Screen by IFA is positive, then CV2 Antibody IgG Titer by IFA will be added. Additional charges apply.

**CPT Codes:** 86255 x2; 84182 ~~x3~~x2; if reflexed, add ~~86256 and/or~~ 84182 x4; ~~if reflexed add~~ 86256 per titer

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

Refer to report

Reference Interval:

Test Number	Components	Reference Interval
	Paraneoplastic Abs (PCCA/ANNA) IgG, CSF	None Detected
	<u>Ma2/Ta Antibody, IgG by Immunoblot, CSF</u>	<u>Negative</u>
	CV2 Ab IgG CBA-IFA Screen, CSF	Less than 1:1
	SOX1 Antibody, IgG by Immunoblot, CSF	Negative
	Amphiphysin Antibody, CSF	Negative

**HOTLINE NOTE:** There is a component change associated with this test. One or more components have been added or removed. Refer to the Hotline Test Mix for interface build information.