

TEST CHANGE

Familial Targeted Sequencing

3005867, FAM NGS

Specimen Requirements:

Patient Preparation:

Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B). New York State Clients: ARUP cannot facilitate testing for New York patients. Please work directly with a New York-approved laboratory.

Specimen Preparation: Transport 2 mL whole blood. (Min: 1 mL)

Transport Temperature: Refrigerated

Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue.

Remarks: Documentation of the familial gene variant from a relative's laboratory test report is required to perform testing. Testing will begin upon receipt of all necessary components, including an original laboratory report detailing the familial variant(s) to be tested.

Stability: Ambient: 72 hours; Refrigerated: 2 weeks; Frozen: Unacceptable

Methodology: Massively Parallel Sequencing

Performed: Varies

Reported: 14-21 days

Note: Documentation of the familial gene variant from a relative's laboratory test report is required to perform testing. Testing will begin upon receipt of all necessary components, including an original laboratory report detailing the familial variant(s) to be tested.

CPT Codes: 81403

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Refer to report

Reference Interval:

By report

HOTLINE NOTE: There is a price change associated with this test. Please contact ARUP Client Services at (800) 522-2787 for additional information.