
New Test **3005867** **Familial Targeted Sequencing** **FAM NGS**



Patient History for Familial Targeted Sequencing Testing



Additional Technical Information

Methodology: Massively Parallel Sequencing
Performed: Varies
Reported: 3 weeks

Specimen Required: Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B).
New York State Clients: ARUP cannot facilitate testing for New York patients. Please work directly with a New York-approved laboratory.
Specimen Preparation: Transport 2 mL whole blood. (Min: 1 mL)
Storage/Transport Temperature: Refrigerated
 Remarks: Documentation of the familial gene variant from a relative's laboratory test report is required to perform testing. Testing will begin upon receipt of all necessary components, including an original laboratory report detailing the familial variant(s) to be tested.
Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue.
Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

Reference Interval: By report

Interpretive Data:
 Refer to report

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA-certified laboratory and is intended for clinical purposes.

Counseling and informed consent are recommended for genetic testing. Consent forms are available online

Note: Documentation of the familial gene variant from a relative's laboratory test report is required to perform testing. Testing will begin upon receipt of all necessary components, including an original laboratory report detailing the familial variant(s) to be tested.

CPT Code(s): 81403

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.