

NEW TEST

[Click for Pricing](#)

Complement C5, Functional
3005960, C5 FUNCT

Specimen Requirements:

Patient Preparation:

Collect: Plain red.

Specimen Preparation: Allow specimen to clot for one hour at room temperature. Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP standard transport tube. (Min: 0.5 mL)

Transport Temperature: CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Separator tubes. Specimens left to clot at 2-8 Degrees C. Specimens exposed to repeated freeze/thaw cycles. Nonfrozen specimens. Grossly hemolyzed or severely lipemic specimens.

Remarks:

Stability: After separation from cells: Ambient: Unacceptable; Refrigerated: 8 days; Frozen: 1 month

Methodology: Quantitative Immunoturbidimetry

Performed: Tue, Fri

Reported: 1-7 days

Note:

CPT Codes: 86161

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

| Component | Interpretation |
|---------------------------|--|
| Complement C5, Functional | Low: 22.9 U/mL or less Low-Normal: 23.0-28.3 U/mL Normal: 28.4 U/mL or greater |

Reference Interval:

| Test Number | Components | Reference Interval |
|-------------|---------------------------|----------------------|
| | Complement C5, Functional | 23.0 U/mL or greater |

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.