

NEW TEST

Click for Pricing

Metabotropic Glutamate Receptor 1 (mGluR1) Antibody, IgG by CBA-IFA With Reflex to Titer, Serum

Effective Date: February 21, 2023

3006044. MGLUR1 SER

3006044, MGLURT SER	
Specimen Requirements:	
Patient Preparation:	
Collect:	Serum separator tube.
Specimen Preparation:	Separate serum from cells within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.2 mL)
Transport Temperature:	Refrigerated.
Unacceptable Conditions:	CSF or plasma. Contaminated, grossly hemolyzed, icteric, or lipemic specimens.
Remarks:	
Stability:	After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month (Three freeze/thaw cycles are
	acceptable)
Methodology:	Semiquantitative Cell-Based Indirect Fluorescent Antibody
Methodology: Performed:	<u> </u>
-	Semiquantitative Cell-Based Indirect Fluorescent Antibody
Performed:	Semiquantitative Cell-Based Indirect Fluorescent Antibody Wed
Performed: Reported:	Semiquantitative Cell-Based Indirect Fluorescent Antibody Wed Within 8 days If mGluR1 antibody IgG is positive, then mGluR1 antibody IgG
Performed: Reported: Note:	Semiquantitative Cell-Based Indirect Fluorescent Antibody Wed Within 8 days If mGluR1 antibody IgG is positive, then mGluR1 antibody IgG titer will be added. Additional charges apply.

Metabotropic glutamate receptor 1 (mGluR1) antibody is found in a subset of patients with autoimmune cerebellar ataxia or autoimmune encephalitis and may occur with or without associated tumor. A negative test result does not rule out a diagnosis of autoimmune cerebellar ataxia or limbic encephalitis. Interpretation of any antineural antibody test requires clinical correlation.

This indirect fluorescent antibody assay utilizes mGluR1 transfected cell lines for detection and



semiquantification of mGluR1 IgG antibody.	
Reference Interval:	
Less than 1:10	

Effective Date: February 21, 2023

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.