

NEW TEST

[Click for Pricing](#)

Hypersensitivity Pneumonitis 1

3006065, HYPER 1

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.15 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma.

Remarks:

Stability: After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)

Methodology: Qualitative Immunodiffusion

Performed: Sun-Sat

Reported: 3-7 days

Note: Testing includes antibodies directed at Aspergillus fumigatus #1, Aspergillus fumigatus #6, Aureobasidium pullulans, pigeon serum, and Micropolyspora faeni.

CPT Codes: 86331 x3; 86606 x2

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval
	A. fumigatus #1 Ab, Precipitin	None detected
	A. fumigatus #6 Ab, Precipitin	None detected
	A. pullulans Ab, Precipitin	None detected
	Pigeon Serum Ab, Precipitin	None detected

	M. faeni Ab, Precipitin	None detected
--	-------------------------	---------------

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.