

NEW TEST

[Click for Pricing](#)

Prenatal Hepatitis Panel

3006343, PRENAT HEP

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST). Also acceptable: Pink (K2EDTA).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 3 mL serum or plasma to an ARUP standard transport tube. (Min: 2.0 mL)

Transport Temperature: Frozen

Unacceptable Conditions: Specimen: Body fluids other than serum or plasma. Condition: Heparinized plasma. Specimens containing particulate material or obvious microbial contamination. Heat-inactivated, severely hemolyzed or lipemic specimens.

Remarks:

Stability: After separation from cells: Ambient: 24 hours; Refrigerated: 5 days; Frozen: 2 months (avoid freeze/thaw cycles)

Methodology: Qualitative Chemiluminescent Immunoassay (CLIA)/Quantitative Transcription-Mediated Amplification (TMA)

Performed: Sun-Sat

Reported: 1-4 days

Note: Order this test only for prenatal specimens. If results for HBsAg screen are reactive (≥ 1.0), then HBsAg Confirmation, Prenatal will be added. Additional charges apply. If the anti-HCV screening result is low positive or high positive, the Hepatitis C Virus by Quantitative NAAT will be added. Additional charges apply.

CPT Codes: 87340; 86803; if reflexed, add 87341; 87522

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

This test should not be used for blood donor screening, associated reentry protocols, or for screening human cells, tissues, and cellular- and tissue-based products (HCT/P).

Components	Interpretation
Hepatitis C Antibody by CIA Index	0.79 IV or less Negative 0.80 to 0.99 IV Equivocal 1.00 to 10.99 IV Low Positive 11.00 IV or greater High Positive

Reference Interval:

Test Number	Components	Reference Interval
	Hepatitis C Antibody by CIA Interp	Negative
	Hepatitis B Surface Antigen, Prenatal	Negative

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.