

NEW TEST

[Click for Pricing](#)

Pharmacogenetics Panel: Psychotropics, with GeneDose Access

3006366, PGXPSYC GD

Specimen Requirements:

Patient Preparation:

Collect: Whole Blood: Lavender (EDTA), pink (K2EDTA), or yellow (ACD solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL)

Transport Temperature: Refrigerated

Unacceptable Conditions: Plasma or serum. Specimens collected in sodium heparin or lithium heparin. Frozen specimens in glass collection tubes.

Remarks:

Stability: Ambient: 72 hours; Refrigerated: 1 week; Frozen: 1 month

Methodology: Polymerase Chain Reaction (PCR)/Fluorescence Monitoring/Sequencing

Performed: Varies

Reported: 5-10 days

Note: Whole blood is the preferred specimen type. Saliva samples that yield inadequate DNA quality and/or quantity will be reported as inconclusive if test performance does not meet laboratory-determined criteria for reporting. Saliva is only validated for the OpenArray and CNV portions of testing and not the long-range PCR/duplication testing. Long-range PCR/duplication testing will not be performed for saliva samples. If long-range PCR/duplication testing is performed, additional charges apply. Approximately less than 5% of samples require 2D6 copy number determination.

CPT Codes: 81225; 81226; 81227; 81230; 81231; 81291; 81479; if reflexed, add 81479

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA-certified laboratory and is intended for clinical purposes.

Counseling and informed consent are recommended for genetic testing. Consent forms are available online.

Reference Interval:

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.