

**NEW TEST – Available Now**

[Click for Pricing](#)

**INSM1 by Immunohistochemistry**

3016567, INSM1 IHC

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Tissue

**Specimen Preparation:** Formalin fix (10 percent neutral buffered formalin) and paraffin embed specimen (cells must be prepared into a cellblock). Protect paraffin block and/or slides from excessive heat. Transport tissue block or 5 unstained (3- to 5-micron thick sections), positively charged slides in a tissue transport kit (recommended but not required). ARUP supply #47808 available online through eSupply using ARUP Connect or contact ARUP Client Services at (800) 522-2787. (Min: 2 slides) If sending precut slides, do not oven bake.

**Transport Temperature:** Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

**Unacceptable Conditions:** Specimens submitted with non-representative tissue type. Depleted specimens.

**Remarks:**

**Stability:** Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

**Methodology:** Immunohistochemistry

**Performed:** Mon-Fri

**Reported:** 1-3 days

**Note:** This test is performed as a stain and return (technical) service only.

**CPT Codes:** 88342

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

Test Number	Components	Reference Interval
-------------	------------	--------------------

**HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.**