

NEW TEST

[Click for Pricing](#)

Cytogenomic SNP Microarray - RAPID

3016596, CMA RAPID

Specimen Requirements:

Patient Preparation: IF Buccal swab: Patient should not eat, drink, smoke or chew gum for 30 minutes before collecting oral sample.

Collect: Peripheral blood or cord blood collected in lavender (K2EDTA). OR Buccal swab: collect one buccal swab using the ORAcollect Collection Kit (ARUP supply #49295). Available online through eSupply using ARUP Connect or contact ARUP Client Services at 800-522-2787

Specimen Preparation: Whole Blood: Transport 3 mL. (Min: 0.5 mL)
Buccal Swab: ensure the sponge tip does not come in contact with any surface prior to collection. Transport buccal swab in ORAcollect Collection Kit.

Transport Temperature: Preferred transport temp: Room temperature.

Unacceptable Conditions: Frozen or clotted specimens

Remarks:

Stability: Whole blood: Room temperature: 2 days; Refrigerated: 3 days;
Frozen: Unacceptable
Buccal swab: Room temperature: 1 week; Refrigerated: 1 week;
Frozen: Unacceptable

Methodology: Genomic Microarray (Oligo-SNP Array)

Note:

CPT Codes: 81229

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Refer to report.

Reference Interval:

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.