

NEW TEST

Click for Pricing

Urine, Amylase

3016646, URN AMY

Specimen Requirements:

Patient Preparation:

Collect:

Specimen Preparation: Urine pH must be adjusted to greater than 7.0. Mix timed urine

collection well. Transfer 4 mL urine to an ARUP standard transport tube. (Min: 0.5 mL) Record total volume, collection time interval, and pH on transport tube or aliquot tube and test

Effective Date: August 21, 2023

request form.

Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens with pH 7.0 or less. Urine collected with acid

Remarks:

Stability: Ambient: 2 days; Refrigerated: 10 days; Frozen: 3 weeks

Methodology: Quantitative Enzymatic Assay

Performed: Sun-Sat

Reported: Within 24 hours

Note:

CPT Codes: 82150

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Random urine: Men: 16-491 U/L Women: 21-447 U/L 24-hour urine: 1 - 17 U/hour

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.