

TEST CHANGE

Platelet Antigen Genotyping Panel, Fetal

3016673, HPAGENO FE

Specimen Requirements:

Patient Preparation:

Collect: Amniotic fluid OR cultured amniocytes OR cultured CVS: Two T-25 flasks at 80 percent confluency. AND ~~maternal whole blood for maternal cell contamination:~~ lavender (~~K2 or K3EDTA~~ EDTA), pink (K2EDTA), or yellow (ACD solution A or B) If the client is unable to culture, order test Cytogenetics Grow and Send (ARUP test code 0040182) in addition to this test and ARUP will culture upon receipt (culturing fees will apply). If you have any questions, contact ARUP's Genetics Processing at 800-522-2787 ext. 3301.

Specimen Preparation: Amniotic fluid: Transport 10 mL amniotic fluid in a sterile container. (Min: 5 mL) OR cultured amniocytes OR cultured CVS: Fill flasks with culture media. Backup cultures must be retained at the client's institution until testing is complete. ~~AND maternal whole~~ Whole blood for maternal cell contamination: Transport ~~2-3~~ mL whole blood. (Min: 1 mL.)

Transport Temperature: Amniotic fluid, cultured amniocytes, cultured CVS: CRITICAL ROOM TEMPERATURE. Must be received within 48 hours of shipment due to lability of cells. ~~Maternal whole~~ Whole blood: ~~room temperature for maternal cell contamination:~~ Refrigerated.

Unacceptable Conditions: Frozen specimens in glass collection tubes.

Remarks: The Patient History Form is available on the ARUP website or by contacting ARUP Client Services.

Stability: Fetal specimens: Ambient: 48 hours; Refrigerated: 48 hours; Frozen: Unacceptable ~~Maternal whole~~ Whole blood ~~for maternal cell contamination:~~ Ambient: 72 hours; Refrigerated: 1 week; Frozen: ~~Unacceptable~~ 1 month

Methodology: Polymerase Chain Reaction (PCR)/Fluorescence Monitoring/Fragment Analysis

Performed: Varies

Reported: 7-14 days

Note:

CPT Codes: 81105; 81106; 81107; 81108; 81109; 81110; 81112; 81265
Fetal Cell Contamination (FCC)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval:
