

NEW TEST

[Click for Pricing](#)

Autoimmune Vision Loss Panel, Serum

3016804, AIVLS

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP standard transport tube. (Min: 0.30 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma. Contaminated, heat-inactivated, hemolyzed, or lipemic specimens.

Remarks:

Stability: After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: Qualitative Immunoblot

Performed: Varies

Reported: 1-8 days

Note: If CV2.1 Antibody IgG Screen by IFA is positive, then CV2.1 Antibody IgG Titer by IFA will be added. Additional charges apply.

CPT Codes: 84182, 86255; if reflexed, add 86256

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Refer to report

Reference Interval:

Test Number	Components	Reference Interval
	CV2.1 Ab IgG CBA-IFA Screen, Serum	Less than 1:10
	Recoverin Ab, IgG by Immunoblot, Serum	Negative

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.