

TEST CHANGE

Autoimmune Vision Loss Panel, Serum

3016804, AIVLS

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP standard transport tube. (Min: 0.30 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma. Contaminated, heat-inactivated, hemolyzed, or lipemic specimens.

Remarks:

Stability: After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: Qualitative Immunoblot/Semi-Quantitative Cell-Based Indirect Fluorescent Antibody

Performed: Varies

Reported: 1-8 days

Note: If CV2-+ Antibody IgG Screen by IFA is positive, then CV2-+ Antibody IgG Titer by IFA will be added. Additional charges apply.

CPT Codes: 84182, 86255; if reflexed, add 86256

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Refer to report

Reference Interval:

Test Number	Components	Reference Interval
	CV2- <u>+</u> Ab IgG CBA-IFA Screen, Serum	Less than 1:10 <u>0</u>
	Recoverin Ab, IgG by Immunoblot, Serum	Negative

HOTLINE NOTE: There is a reflexive pattern change associated with this test. One or more orderable or component has been added or removed to the reflexive pattern. Refer to the Hotline Test Mix for interface build information.