

NEW TEST - Available Now

Click for Pricing

Adiponectin, Quantitative Serum/Plasma

3017195, ADIP SP

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST), light green (lithium heparin with

gel separator), or plain red.

Specimen Preparation: Allow serum to clot for 15-20 minutes at room temperature.

Separate serum or plasma from cells ASAP or within 2 hours of

Effective Date: May 3, 2024

collection. Transfer 0.5 mL serum or plasma to an ARUP

standard transport tube. (Min: 0.2 mL)

Transport Temperature: Frozen

Unacceptable Conditions: EDTA plasma.

Remarks:

Stability: After separation from cells: Ambient: 2 days; Refrigerated: 1

week; Frozen: 1 month

Methodology: Quantitative Radioimmunoassay (RIA)

Performed: Wed

Reported: 1-8 days

Note:

CPT Codes: 83519

New York DOH Approval Status: Specimens from New York clients will be sent out to a New

York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

Effective Date: May 3, 2024



Age Male (ug/mL) Female(ug/mL) 0-7 years 2.33-26.5 2.33-26.5 8-9 years 3.96-14.9 3.96-14.9 10-11 years 3.36-13.8 3.36-13.8 12-13 years 4.50-13.2 4.50-13.2 14-15 years 3.67-13.7 3.67-13.7 16-19 years 2.74-13.3 2.74-13.3 greater than 19 2.00-13.9 4.00-19.4 years

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.