

NEW TEST

[Click for Pricing](#)

Human Immunodeficiency Virus 1 and 2 (HIV-1/HIV-2) by Qualitative NAAT

3017779, HIV1,2QUAL

Specimen Requirements:

Patient Preparation:

Collect: Lavender (EDTA), pink (K2EDTA), plasma preparation tube (PPT), or serum separator tube (SST).

Specimen Preparation: Separate from cells within 24 hours of collection. Transfer 2 mL plasma to an ARUP standard transport tube (ARUP supply #15824). Available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787. (Minimum volume: 1.3mL)

Transport Temperature: Frozen

Unacceptable Conditions: Heparinized specimens.

Remarks:

Stability: After separation from cells: Ambient: 24 hrs; Refrigerated: 5 days; Frozen: 6 weeks.

Methodology: Qualitative Polymerase Chain Reaction (PCR)

Performed: Sun-Sat

Reported: 1-4 days

Note: Test detects and differentiates HIV-1 and HIV-2 virus RNA. Proviral DNA will not be detected.

CPT Codes: 87535

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.



*A nonprofit enterprise of the University of Utah
and its Department of Pathology*

Effective Date: **October 21, 2024**