

NEW TEST

[Click for Pricing](#)

Gabapentin, Serum or Plasma

3017893, GABAP-SP

Specimen Requirements:

Patient Preparation:

Collect: Plain red. Also acceptable: Dark green (sodium or lithium heparin), lavender (K2 or K3EDTA) or pink (K2EDTA).

Specimen Preparation: Separate serum or plasma from cells within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP standard transport tube. (Min: 0.2 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Whole blood. Gel separator tubes, light blue (citrate), or yellow (SPS or ACD solution).

Remarks:

Stability: After separation from cells: Ambient: 48 Hours; Refrigerated: 1 week; Frozen: 1 month 3 freeze/thaw cycles

Methodology: Quantitative Enzyme Immunoassay (EIA)

Performed: Sun-Sat

Reported: 1-2 days

Note:

CPT Codes: 80171

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.