

NEW TEST – Available Now

[Click for Pricing](#)

Vasoactive Intestinal Polypeptide (VIP), Plasma

3018202, VIP PLA

Specimen Requirements:

Patient Preparation:

Collect: Chilled lavender or pink (K2EDTA or K3EDTA).

Specimen Preparation: Collect in a prechilled tube. Mix well, then place on ice until centrifugation. Transfer 2 mL of plasma to an ARUP standard transport tube (Min: 0.5 mL). Freeze at -20.

Transport Temperature: Frozen. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Grossly hemolyzed, lipemic, icteric, or clotted specimens.

Remarks:

Stability: After separation from cells: Room Temperature: 4 hours; Refrigerated: 24 hours; Frozen: 3 months.

Methodology: Quantitative Enzyme-Linked Immunosorbent Assay (ELISA)

Performed: Mon, Thu

Reported: 3-7 days

Note:

CPT Codes: 84586

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

Test Number	Test Name	Reference Interval
3018312	Vasoactive Intestinal Polypeptide (VIP)	0-89.1 pg/mL

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.



*A nonprofit enterprise of the University of Utah
and its Department of Pathology*

Effective Date: **October 21, 2024**