

## **NEW TEST - Available Now**

Click for Pricing

Kelch-Like Protein 11 Antibody, IgG by CBA-IFA, With Reflex to Titer, CSF 3018508, KLHL11 CSF

Effective Date: January 21, 2025

Specimen Requirements:	
Patient Preparation:	
Collect:	Separate CSF.
Specimen Preparation:	Transfer 0.5 mL CSF to an ARUP standard transport tube. (Min: 0.15 mL)
Transport Temperature:	Refrigerated
Unacceptable Conditions:	Grossly hemolyzed or contaminated specimens.
Remarks:	
Stability:	Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month (three freeze/thaw cycles are acceptable)
	(
Methodology:	Semi-Quantitative Cell-Based Indirect Fluorescent Antibody
Methodology: Performed:	
J.	Semi-Quantitative Cell-Based Indirect Fluorescent Antibody
Performed:	Semi-Quantitative Cell-Based Indirect Fluorescent Antibody Wed
Performed: Reported:	Semi-Quantitative Cell-Based Indirect Fluorescent Antibody  Wed  1-8 days  If KLHL11 antibody IgG is positive, then KLHL11 antibody IgG
Performed: Reported: Note:	Semi-Quantitative Cell-Based Indirect Fluorescent Antibody  Wed  1-8 days  If KLHL11 antibody IgG is positive, then KLHL11 antibody IgG titer will be added. Additional charges apply.

Interpretive Data:

IgG antibodies to KLHL11 are associated with paraneoplastic neurologic syndromes with phenotypes most often including a combination of brainstem and cerebellar encephalitis as well as sensorineural hearing loss. Patients with anti-KLHL11 syndromes should be thoroughly evaluated for cancer, including testicular cancer, as neurologic symptoms often precede cancer diagnosis. Consider sending testing in serum as well as CSF to improve diagnostic yield. Coexisting and clinically relevant antineural antibodies have been reported; consider ordering a phenotype-specific panel to assess for these. Results (positive or negative) should be interpreted in the context of the patient's complete clinical picture, as false positives may occur, and a negative result does not exclude the diagnosis of immune-mediated neurologic disease.

Reference Interval:



Test Components Reference Interval

KLHL11 Ab IgG CBA-IFA Screen, CSF Less than 1:1

Effective Date: January 21, 2025

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.