

NEW TEST

[Click for Pricing](#)

Tobramycin Level, Random, Serum

3018760, TOBRAM R

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST). Also acceptable: Plain red

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 0.5 mL serum to an ARUP standard transport tube. (Min: 0.25 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen. Also acceptable: Refrigerated

Unacceptable Conditions:

Remarks:

Stability: Ambient: 72 hours; Refrigerated: 1 week; Frozen: 28 days

Methodology: Quantitative Enzyme Immunoassay (EIA)

Performed: Varies

Reported: 3-5 days

Note:

CPT Codes: 80200

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.