

NEW TEST

[Click for Pricing](#)

Polymyositis Panel

3018868, POLY MY02

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST)

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer two 1 mL serum aliquots to ARUP standard transport tubes. (Min: 0.5 mL/aliquot)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Hemolyzed, hyperlipemic, icteric, heat-treated, or contaminated specimens

Remarks:

Stability: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: Qualitative Immunoprecipitation / Semi-Quantitative Multiplex Bead Assay / Semi-Quantitative Indirect Fluorescent Antibody (IFA) / Qualitative Immunoblot

Performed: Mon, Tue, Thu, Fri

Reported: 7-18 days

Note: Antibodies: PL-7, PL12, EJ, OJ, SRP, Jo-1, ANA, Ha, Ks, Zo

CPT Codes: 83516 x5; 86235; 84182 x3; 86039

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.