

**TEST CHANGE**

**Interstitial Lung Disease Autoantibody Panel**

3018869, ILD PANEL2

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Serum separator tube (SST).

**Specimen Preparation:** Separate from cells ASAP or within 2 hours of collection. Transfer ~~3 to 4~~<sup>five</sup> ~~1~~ mL ~~serum aliquots~~ to ARUP standard transport tubes. (Min: ~~1.5~~<sup>2</sup>~~8~~ mL)

**Transport Temperature:** Refrigerated.

**Unacceptable Conditions:** Hemolyzed, hyperlipemic, icteric, heat-treated, or contaminated specimens.

**Remarks:**

**Stability:** Ambient: 24 hours; Refrigerated: 1 weeks; Frozen: 1 month

**Methodology:** Qualitative Immunoprecipitation / Semi-Quantitative Multiplex Bead Assay / Qualitative Immunoblot / Semi-Quantitative Enzyme-Linked Immunosorbent Assay (ELISA) / Quantitative Immunoturbidimetry / Semi-Quantitative Indirect Fluorescent Antibody (IFA)  
[/ Qualitative Particle-Based Multianalyte Technology \(PMAT\)](#)

**Performed:** Sun-Sat

**Reported:** 7-18 days

**Note:** Antibodies: Ro52, Ro60, Jo-1, PL-7, PL12, EJ, Ku, SRP, OJ, PM/Scl-100, MDA5, CCP, Scl-70, RA, ANA, NXP-2, RNA Polymerase III, Ha, Ks, Zo

**CPT Codes:** 86235 x5; 83516 ~~x9~~<sup>x7</sup>; 84182 ~~x3~~<sup>x5</sup>; 86431; 86200; 86039

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

SSA-52 (Ro52) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive
SSA-60 (Ro60) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive
Scleroderma (Scl-70) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive
Jo-1 Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive
Cyclic Citrullinated Peptide (CCP) Ab, IgG/A	19 Units or less: Negative 20-39 Units: Weak positive 40-59 Units: Moderate positive 60 Units or Greater: Strong positive
RNA Polymerase III Antibody, IgG	19 Units or less: Negative 20-39 Units: Weak positive 40-80 Units: Moderate positive 81 Units or Greater: Strong positive

Reference Interval:

Test Number	Components	Reference Interval
	Zo (phenylalanyl-tRNA synthetase) Ab	Negative
	Ks (asparaginyl-tRNA synthetase) Ab	Negative
	Ha (tyrosyl-tRNA synthetase) Ab	Negative
	Cyclic Citrullinated Peptide Ab, IgG/A	19 Units or less
	NXP2 (Nuclear matrix protein-2) Ab	Negative
	MDA5 (CADM-140) Ab	Negative
	Antinuclear Antibody (ANA), HEp-2, IgG	Less than 1:80
	SSA-60 (Ro60) (ENA) Antibody, IgG	40 AU/mL or less
	OJ (isoleucyl-tRNA synthetase) Antibody	Negative
	SRP (Signal Recognition Particle) Ab	Negative
	Ku Antibody	Negative
	EJ (glycyl-tRNA synthetase) Antibody	Negative
	PL-12 (alanyl-tRNA synthetase) Antibody	Negative
	PL-7 (threonyl-tRNA synthetase) Antibody	Negative
	PM/Scl 100 Antibody, IgG	Negative
	RNA Polymerase III Antibody, IgG	19 Units or less
	Jo-1 (Histidyl-tRNA Synthetase) Ab, IgG	40 AU/mL or less
	SSA-52 (Ro52) (ENA) Antibody, IgG	40 AU/mL or less
	Scleroderma (Scl-70) (ENA) Antibody, IgG	40 AU/mL or less
	Rheumatoid Factor	0-14 IU/mL