

**NEW TEST – Available Now**

[Click for Pricing](#)

**HBV Prenatal Triple Panel**

3018943, HBV TRI PN

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Serum separator tube (SST).

**Specimen Preparation:** Separate serum from cells ASAP or within 2 hours of collection. Transfer 3.00 mL serum to an ARUP standard transport tube. (Min: 1.750 mL)

**Transport Temperature:** Refrigerated.

**Unacceptable Conditions:** Heparinized plasma. Specimens containing particulate material or obvious microbial contamination. Heat-inactivated or severely hemolyzed.

**Remarks:**

**Stability:** After separation from cells: Ambient: 12 hours; Refrigerated: 7 days; Frozen: 30 days (avoid repeated freeze/thaw cycles).

**Methodology:** Quantitative Chemiluminescent Immunoassay (CLIA) / Qualitative Chemiluminescent Immunoassay (CLIA)

**Performed:** Sun-Sat

**Reported:** 1-2 days

**Note:** The HBV core total assay tests for IgG and IgM antibodies, but does not differentiate between them. HBsAb results greater than 1,000.00 IU/L are reported as greater than 1,000.00 IU/L. If results for HBsAg screen are reactive (=1.0), then HBsAg Confirmation, Prenatal will be added. Additional charges apply.

**CPT Codes:** 86706; 86704; 87340; if reflexed, add 87341

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

Test Number	Components	Reference Interval	
	Hepatitis B Surface Antigen, Prenatal	Negative	
	Hepatitis B Surface Antigen, Prenatal		
	Hepatitis B Surface Antigen, Prenatal		
	Hepatitis B Core Antibodies, Total	Negative	
	Hepatitis B Surface Antibody		
		Less than 10.00 IU/L	Negative
		Greater than or equal to 10.00 IU/L	Positive

**HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.**