

## NEW TEST – Available Now

[Click for Pricing](#)

### Measles Virus by Qualitative NAAT

3019269, MEASLESPCR

#### Specimen Requirements:

##### Patient Preparation:

**Collect:** Preferred: Respiratory swab (nasal, nasopharyngeal, oropharyngeal, and throat) Also Acceptable: Urine in Aptima Urine Collection Kit.

**Specimen Preparation:** Swab: Place swab in viral transport media (ARUP Supply #12884). Available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787. Urine: Transfer 2 mL urine within 24 hours to an Aptima Urine Specimen Transport Tube (ARUP supply #28908 PK/50 or 54556 PK/10). Liquid level must be between fill lines on tube.

**Transport Temperature:** Frozen

**Unacceptable Conditions:** Swabs not in viral transport media.

##### Remarks:

**Stability:** Ambient: 3 days; Refrigerated: 7 days; Frozen: 30 days

**Methodology:** Qualitative Nucleic Acid Amplification Test (NAAT)

**Performed:** Sun-Sat

**Reported:** 1-3 days

##### Note:

**CPT Codes:** 87798 x2

**New York DOH Approval Status:** Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

##### Interpretive Data:

A negative result does not rule out the presence of PCR inhibitors in the patient specimen or assay-specific nucleic acid in concentrations below the level of detection by this assay.

##### Reference Interval:

Test Number	Components	Reference Interval
	Measles Virus, NAAT	Not Detected
	Measles Virus Vaccine Strain, NAAT	Not Detected

**HOTLINE NOTE:** Refer to the Hotline Test Mix for interface build information.