

NEW TEST – Available Now

[Click for Pricing](#)

Alpha Globin (HBA1 and HBA2) Sequencing and Deletion/Duplication, Fetal

3019566, HBA FGA FE

Specimen Requirements:

Patient Preparation:

Collect: Fetal Specimen: Two T-25 flasks at 90% confluent of cultured amniocytes or cultured chorionic villus sampling (CVS). AND Maternal Whole Blood Specimen: Lavender (EDTA), pink (K2EDTA), or yellow (ACD solution A or B).

Specimen Preparation: Cultured Amniocytes or Cultured CVS: Fill flasks with culture media. Transport two T-25 flasks at 90 percent confluent of cultured amniocytes or cultured CVS filled with culture media. Backup cultures must be retained at the client's institution until testing is complete. If ARUP receives a sample below the minimum confluence, Cytogenetics Grow and Send (0040182) will be added on by ARUP, and additional charges will apply. If clients are unable to culture specimens, Cytogenetics Grow and Send should be added to initial order. Maternal Whole Blood Specimen: Transport 3 mL whole blood. (Min: 1 mL).

Transport Temperature: Cultured Amniocytes or Culture CVS: CRITICAL ROOM TEMPERATURE. Separate specimens must be submitted when multiple tests are ordered. Must be received within 48 hours of shipment due to viability of cells. Maternal Specimen: Refrigerated

Unacceptable Conditions: Frozen specimens

Remarks:

Stability: Cultured Amniocytes or Cultured CVS: Room temperature: 48 hours; Refrigerated: Unacceptable; Frozen: Unacceptable
Maternal Whole Blood Specimen: Room temperature: 7 days; Refrigerated: 1 month; Frozen: Unacceptable

Methodology: Qualitative Multiplex Ligation-Dependent Probe Amplification (MLPA)/Qualitative Sequencing

Performed: Varies

Reported: 14-21 days

Note: Reported times are based on receiving the two T-25 flasks at 90 percent confluent. Cell culture time is independent of

testing turnaround time. Maternal specimen is recommended for proper test interpretation. Order Maternal Cell Contamination.

CPT Codes: 81259; 81269; 81265

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report

Reference Interval:

By report

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.