

NEW TEST

[Click for Pricing](#)

DNA Extract and Hold for Whole Blood and Bone Marrow

3020699, WBBMDNAEXT

Specimen Requirements:

Patient Preparation:

Collect: Whole blood or bone marrow in lavender (K3EDTA or K2EDTA) or pink (K2EDTA). Also Acceptable: Whole blood or bone marrow in green (sodium heparin) or yellow (ACD).

Specimen Preparation: Do not freeze. Transport 1 mL whole blood or bone marrow.

Transport Temperature: Preferred transport temp: Refrigerated

Unacceptable Conditions: Unacceptable specimens types: Plasma, serum, FFPE tissue blocks/slides, or fresh or frozen tissue, bone core, white cell pellet, other body fluids.
Unacceptable collected device: Any anticoagulants other than EDTA, sodium heparin, or yellow.
Unacceptable specimen conditions: Clotted or grossly hemolyzed specimens.

Remarks:

Stability: Room Temperature: 3 days; Refrigerated: 30 days; Frozen: Unacceptable

Methodology: DNA Extraction

Note: Yield of 50 microliters. DNA will be held for 3 months for possible add-on testing. Additional charges will apply if extracted DNA is sent to an outside laboratory.

CPT Codes: NA

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.