

Apolipoprotein E Genotyping, Cardiovascular Disease Risk

The *APOE* gene encodes the apolipoprotein E (apoE) protein, which combines with lipids to form lipoproteins that transport cholesterol and other fats through the bloodstream.

The three most common *APOE* variants are designated as the e2, e3, and e4 alleles. E3 is the wild type allele, e4 is associated with [Alzheimer disease](#), and e2 is associated with an increased risk for early cardiovascular disease (CVD) and hyperlipoproteinemia type III (HLP III).

HLP III accounts for up to 5% of premature coronary heart disease (CHD) and may also present with elevated cholesterol, triglycerides, and very low density lipoprotein (VLDL). Early identification allows for treatment with lipid-lowering agents.

Genetics

Gene

APOE

Inheritance of Hyperlipoproteinemia Type III

Multifactorial

Penetrance

1-5% percent of individuals homozygous for the e2 allele develop HPL III

Structure

ApoE is a critical protein component of VLDL and chylomicrons

Variants

Three common alleles (e2, e3, e4) differ at amino acid positions 112 (130 legacy) and 158 (176 legacy)

- Allele frequencies:
 - e2 (c.388T; p.130Cys and c.526C>T; p.Arg176Cys): 10%
 - e3 (c.388T; p.130Cys and c.526C; p.176Arg): 75%
 - e4 (c.388T>C; p.Cys130Arg and c.526C; p.176Arg): 15%
- *APOE* e2: binds the lipoprotein receptors with only 2% of the affinity of e3 and e4 isoforms
 - Results in impaired clearance of chylomicron and VLDL remnants
 - Leads to increased plasma cholesterol and triglyceride levels
- Homozygosity for e2 is present in 1% of White individuals
 - Only genotype associated with HLP III
 - 1-4% of homozygotes will develop HLP III
 - Found in >90% of individuals with HLP III
- *APOE* e3: considered wild type
- *APOE* e4: associated with increased plasma cholesterol

Test Interpretation

Sensitivity/Specificity

- Clinical sensitivity: >90% for individuals with HLP III¹
- Analytical sensitivity/specificity: 99%

Featured ARUP Testing

[Apolipoprotein E \(APOE\) Genotyping, Cardiovascular Risk 2013337](#)

Method: Polymerase Chain Reaction/Fluorescence Monitoring

- Determines *APOE* genotype in context of HLP III for evaluation of premature CHD
- Use for cardiovascular risk assessment only
- Not recommended for asymptomatic individuals <18 years

Results

- *APOE* e2/e2: provides additional evidence for a clinical diagnosis of HLP III; by itself, genotype is not diagnostic for HLP III
- *APOE* e3/e3: most common genotype found in general population
- *APOE* e4/e4: associated with increased plasma cholesterol levels that may contribute to CHD
- *APOE* e2/e3, e2/e4, e3/e4: no significantly increased risk for HLP III
- *APOE* e2/e4 and e3/e4: some association with increased plasma cholesterol levels and atherosclerosis

Limitations

- Diagnostic errors can occur due to rare sequence variations
- Rare *APOE* variants and variants in other genes that cause HLP III are not detected
- *APOE* e2 homozygosity is neither sufficient nor necessary to cause HPL III

References

1. Eichner JE, Dunn T, Perveen G, et al. [Apolipoprotein E polymorphism and cardiovascular disease: a HuGE review](#). *Am J Epidemiol*. 2002;155(6):487-495.

Related Information

[Apolipoprotein E \(APOE\) Genotyping, Alzheimer Disease Risk](#)

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology. 500 Chipeta Way, Salt Lake City, UT 84108
(800) 522-2787 | (801) 583-2787 | [aruplab.com](#) | [arupconsult.com](#)
Content Review September 2020 | Last Update June 2023