

Chronic Lymphocytic Leukemia Mutation Panel by Next Generation Sequencing

Chronic lymphocytic leukemia (CLL) is a hematopoietic disorder characterized by chronic monoclonal B-cell proliferation. Recent studies have identified recurrently mutated genes with diagnostic and/or prognostic impact in CLL and other lymphoid malignancies. The presence of certain mutations may inform clinical management. This multigene panel by massively parallel sequencing (next generation sequencing) is a more cost-effective approach when compared to the cost of multiple single gene tests. This test can be used to complement the morphologic and cytogenetic workup of CLL and other lymphoid malignancies.

Disease Overview

Diagnostic Issues

- Genetic targets contained in the panel are relevant in CLL and other lymphoid malignancies
- Identification of one or more clonal genetic abnormalities may aid in establishing the diagnosis of a neoplasm, but results must be interpreted within the context of other clinical and hematologic findings
- Identification of certain variants or patterns of variants may aid in prognostication of patients with CLL

Prognostic and Treatment Issues

- Certain variants or patterns of variants may have prognostic significance
- Certain variants may be sensitive to or provide resistance to targeted therapies

Genetics

Genes Tested

ATM, BCL2, BIRC3^a, BRAF, BTG1, BTK, CARD11, CD79B, CXCR4, DDX3X, FBXW7, IKZF3, KRAS, MAP2K1, MED12, MGA, MYD88, NOTCH1, NRAS, PLCG2, POT1, RPS15^a, SAMHD1, SF3B1, TP53, XPO1, ZMYM3

(^aOne exon of the preferred transcript is not covered by sequencing for the indicated gene. See [Genes Tested](#) table below for full list of targeted regions and exclusions.)

Test Interpretation

Results

- Positive: a somatic variant in one of the tested genes was detected
 - Clinical relevance will be described, if known
- Negative: no variants were detected in the sequenced genes

Tests to Consider

[Chronic Lymphocytic Leukemia Mutation Panel by Next Generation Sequencing 3001858](#)

Method: Massively Parallel Sequencing

See [Related Tests](#)

For more information on individual related tests, see [Chronic Lymphocytic Leukemia](#).

Limitations

- A negative result does not exclude a diagnosis of CLL
- Not intended to detect minimal residual disease (MRD)
- This test does not determine *IGHV* mutation status
- Variants may be present below the limit of detection (LOD) of 5% allele frequency
- Variants greater than 24 base pairs may be detected at LOD, but the analytical sensitivity may be reduced
- Variants may not be identified due to technical limitations in the presence of pseudogenes or in repetitive or homologous regions
- Variants in regions that are not included in the preferred transcript for the targeted genes will not be detected; see [Genes Tested](#) table below for full list of targeted regions and exclusions
- Interpretation of this test result may be impacted if this patient has had an undisclosed allogenic bone marrow transplant or stem cell transplant
- Does not detect translocations, gene rearrangements, copy number alterations, or microsatellite instability
- Does not distinguish between somatic and germline variants

Analytical Sensitivity

Variant Class	Analytical Sensitivity (PPA) ^a Estimate (%)	Analytical Sensitivity (PPA) 95% Credibility Region ^a (%)
SNVs	96.9%	95.1-98.1%
Insertions/Duplications (1-24bp)	98.1%	95.5-99.3%
Insertions/Duplications (>24bp)	> 99%	92.9-100.0%
Deletions (1-24bp)	96.7%	92.8-98.7%
Deletions (>24bp)	90%	79.5-96.1%
MNVs	97%	93.0-99.0%

^aGenes included on this test are a subset of a larger methods-based validation from which the PPA values are derived.

bp, base pairs; MNVs, multinucleotide variants; PPA, positive percent agreement; SNVs, single nucleotide variants

Genes Tested	
Gene	Preferred Transcript ^{a, b}
<i>ATM</i>	NM_000051
<i>BCL2</i>	NM_000633
<i>BIRC3</i>	NM_001165 (exon 5 excluded)
<i>BRAF</i>	NM_004333

^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the NOTCH1 3'UTR. In addition, coding exons noted here as excluded are not sequenced due to technical limitations of the assay.

Gene	Preferred Transcript ^{a, b}
<i>BTG1</i>	NM_001731
<i>BTK</i>	NM_000061
<i>CARD11</i>	NM_032415
<i>CD79B</i>	NM_000626
<i>CXCR4</i>	NM_003467
<i>DDX3X</i>	NM_001193416
<i>FBXW7</i>	NM_033632
<i>IKZF3</i>	NM_012481
<i>KRAS</i>	NM_004985
<i>MAP2K1</i>	NM_002755
<i>MED12</i>	NM_005120
<i>MGA</i>	NM_001164273
<i>MYD88</i>	NM_002468
<i>NOTCH1</i>	NM_017617
<i>NRAS</i>	NM_002524
<i>PLCG2</i>	NM_002661
<i>POT1</i>	NM_015450
<i>RPS15</i>	NM_001018 (exon 19 excluded)
<i>SAMHD1</i>	NM_015474
<i>SF3B1</i>	NM_012433
<i>TP53</i>	NM_000546
<i>XPO1</i>	NM_003400
<i>ZMYM3</i>	NM_201599

^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the NOTCH1 3'UTR. In addition, coding exons noted here as excluded are not sequenced due to technical limitations of the assay.

Related Information

[Chronic Lymphocytic Leukemia - CLL](#)

Related Tests

[Chromosome FISH, CLL Panel 2002295](#)

Method: Fluorescence in situ Hybridization

[IGHV Mutation Analysis by Sequencing 0040227](#)

Method: Polymerase Chain Reaction/Sequencing

[Leukemia/Lymphoma Phenotyping Evaluation by Flow Cytometry 3001780](#)

Method: Flow Cytometry

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology. 500 Chipeta Way, Salt Lake City, UT 84108
| (800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com
Content Review December 2019 | Last Update March 2020