

Myeloid Malignancies and Acute Myeloid Leukemia Mutation Panels by Next Generation Sequencing

Myeloid malignancies are clonal disorders of hematopoietic stem and progenitor cells that include myelodysplastic syndromes (MDS), myeloproliferative neoplasms (MPN), myelodysplastic/myeloproliferative neoplasms (MDS/MPN), and acute myeloid leukemia (AML). Recent studies have identified recurrently mutated genes with diagnostic and/or prognostic impact in myeloid malignancies. The presence of certain mutations may inform clinical management. This multigene panel by massively parallel sequencing (next generation sequencing) is a more cost-effective approach when compared to the cost of multiple single gene tests. This test can be used to complement the morphologic and cytogenetic workup of myeloid malignancies.

Disease Overview

Diagnostic Issues

- Genetic targets contained in panels are relevant across the spectrum of myeloid malignancies
- Identification of one or more clonal genetic abnormalities may aid in establishing the diagnosis and subclassification of a myeloid neoplasm
- Identification of certain variants or patterns of variants may aid in prognostication and clinical management of patients with a diagnosis of myeloid malignancy

Prognostic and Treatment Issues

- Certain variants or patterns of variants may have diagnostic or prognostic significance
- Certain variants may inform clinical management

Genetics

Genes Tested: Myeloid Malignancies Mutation Panel by Next Generation Sequencing

ANKRD26, ASXL1, ASXL2, BCOR, BCORL1, BRAF, CALR, CBL, CBLB, CEBPA, CSF3R, CUX1^a, DDX41, DNMT1^a, DNMT3A, ELANE, ETNK1, ETV6, EZH2, FBXW7, FLT3, GATA1, GATA2, GNAS, HNRNP, IDH1, IDH2, IL7R, JAK1, JAK2, JAK3, KDM6A^a, KIT, KMT2A, KRAS, LUC7L2, MPL, NOTCH1, NPM1^a, NRAS, NSD1, PHF6, PIGA, PRPF40B, PRPF8, PTPN11, RAD21, RUNX1, SETBP1, SF3B1, SH2B3, SMC1A, SMC3, SRSF2, STAG2, STAT3, STAT5B^a, SUZ12^a, TET2, TP53, U2AF1, U2AF2, WT1, ZRSR2

Genes Tested: Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing

ANKRD26, ASXL1, CEBPA, DDX41, DNMT3A, ETV6, FLT3, GATA2, IDH1, IDH2, KIT, KRAS, NPM1^a, NRAS, RUNX1, TP53, WT1

(^aOne or more exons of the preferred transcript are not covered by sequencing for the indicated gene. See [Genes Tested](#) table below for full list of targeted regions and exclusions.)

Test Interpretation

Tests to Consider

[Myeloid Malignancies Mutation Panel by Next Generation Sequencing 2011117](#)

Method: Massively Parallel Sequencing

[Myeloid Malignancies Somatic Mutation and Copy Number Analysis Panel 2012182](#)

Method: Massively Parallel Sequencing/Genomic Microarray (Oligo-SNP Array)

[Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing 3002714](#)

Method: Massively Parallel Sequencing

For more information on genomic microarray testing in oncology, see the additional technical information document, [Cytogenomic Microarray – Oncology](#)

See Related Tests

For additional test information, including information on individual tests, refer to the [Acute Myeloid Leukemia Molecular Genetic Testing](#) Test Fact Sheet.

Results



- Positive: a somatic variant in one of the tested genes was detected
 - Clinical relevance will be described, if known
- Negative: no variants were detected in the sequenced genes

Limitations

- Not intended to detect minimal residual disease (MRD)
- Variants may be present below the limit of detection (LOD) of 5% allele frequency
- Variants greater than 24 base pairs may be detected at LOD, but the analytical sensitivity may be reduced
- Variants may not be identified due to technical limitations in the presence of pseudogenes or in repetitive or homologous regions
- Variants in regions that are not included in the preferred transcript for the targeted genes will not be detected; see [Genes Tested](#) table below for full list of targeted regions and exclusions

Analytical Sensitivity

Variant Class	Analytical Sensitivity (PPA) ^a Estimate (%)	Analytical Sensitivity (PPA) 95% Credibility Region ^a (%)
SNVs	96.9%	95.1-98.1%
Insertions/Duplications (1-24bp)	98.1%	95.5-99.3%
Insertions/Duplications (>24bp)	> 99%	92.9-100.0%
Deletions (1-24bp)	96.7%	92.8-98.7%
Deletions (>24bp)	90%	79.5-96.1%
MNVs	97%	93.0-99.0%
FLT3 ITDs	>99%	97.1-100.0%

^aGenes included on this test are a subset of a larger methods-based validation from which the PPA values are derived.

bp, base pairs; ITDs, internal tandem duplications; MNVs, multinucleotide variants; PPA, positive percent agreement; SNVs, single nucleotide variants

Genes Tested by Myeloid Malignancies Mutation Panel by Next Generation Sequencing

(The subset of genes tested by Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing are bolded)

Gene	Preferred Transcript ^a	Excluded Exons ^b
ANKRD26	NM_014915	
ASXL1	NM_015338	
ASXL2	NM_018263	
BCOR	NM_001123385	
BCORL1	NM_021946	
BRAF	NM_004333	
CALR	NM_004343	



^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the ANKRD26 5'UTR and NOTCH1 3'UTR. In addition, coding exons noted here are not sequenced due to technical limitations of the assay.

Gene	Preferred Transcript ^a	Excluded Exons ^b
<i>CBL</i>	NM_005188	
<i>CBLB</i>	NM_170662	
<i>CEBPA</i>	NM_004364	
<i>CSF3R</i>	NM_156039	
<i>CUX1</i>	NM_181552	24
<i>DDX41</i>	NM_016222	
<i>DNMT1</i>	NM_001130823	5
<i>DNMT3A</i>	NM_175629	
<i>ELANE</i>	NM_001972	
<i>ETNK1</i>	NM_018638	
<i>ETV6</i>	NM_001987	
<i>EZH2</i>	NM_004456	
<i>FBXW7</i>	NM_033632	
<i>FLT3</i>	NM_004119	
<i>GATA1</i>	NM_002049	
<i>GATA2</i>	NM_032638	
<i>GNAS</i>	NM_000516	
<i>HNRNPK</i>	NM_002140	
<i>IDH1</i>	NM_005896	
<i>IDH2</i>	NM_002168	
<i>IL7R</i>	NM_002185	
<i>JAK1</i>	NM_002227	
<i>JAK2</i>	NM_004972	
<i>JAK3</i>	NM_000215	
<i>KDM6A</i>	NM_001291415	13



^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the ANKRD26 5'UTR and NOTCH1 3'UTR. In addition, coding exons noted here are not sequenced due to technical limitations of the assay.

Gene	Preferred Transcript ^a	Excluded Exons ^b
<i>KIT</i>	NM_000222	
<i>KMT2A</i>	NM_001197104	
<i>KRAS</i>	NM_004985	
<i>LUC7L2</i>	NM_016019	
<i>MPL</i>	NM_005373	
<i>NOTCH1</i>	NM_017617	
<i>NPM1</i>	NM_002520	1
<i>NRAS</i>	NM_002524	
<i>NSD1</i>	NM_022455	
<i>PHF6</i>	NM_001015877	
<i>PIGA</i>	NM_002641	
<i>PRPF40B</i>	NM_001031698	
<i>PRPF8</i>	NM_006445	
<i>PTPN11</i>	NM_002834	
<i>RAD21</i>	NM_006265	
<i>RUNX1</i>	NM_001754	
<i>SETBP1</i>	NM_015559	
<i>SF3B1</i>	NM_012433	
<i>SH2B3</i>	NM_005475	
<i>SMC1A</i>	NM_006306	
<i>SMC3</i>	NM_005445	
<i>SRSF2</i>	NM_003016	
<i>STAG2</i>	NM_001042749	
<i>STAT3</i>	NM_139276	
<i>STAT5B</i>	NM_012448	6-9



^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the ANKRD26 5'UTR and NOTCH1 3'UTR. In addition, coding exons noted here are not sequenced due to technical limitations of the assay.

Gene	Preferred Transcript ^a	Excluded Exons ^b
<i>SUZ12</i>	NM_015355	1-9
<i>TET2</i>	NM_001127208	
<i>TP53</i>	NM_000546	
<i>U2AF1</i>	NM_006758	
<i>U2AF2</i>	NM_007279	
<i>WT1</i>	NM_024426	
<i>ZRSR2</i>	NM_005089	

^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the ANKRD26 5'UTR and NOTCH1 3'UTR. In addition, coding exons noted here are not sequenced due to technical limitations of the assay.

Related Information

[Acute Myeloid Leukemia - AML](#)
[Myelodysplastic Syndromes](#)
[Myeloproliferative Neoplasms](#)

Related Tests

[CEBPA Mutation Detection 2004247](#)

Method: Polymerase Chain Reaction/Sequencing

[NPM1 Mutation Detection by RT-PCR, Quantitative 3000066](#)

Method: Quantitative Reverse-Transcription Polymerase Chain Reaction

[IDH1 and IDH2 Mutation Analysis, exon 4 2006444](#)

Method: Polymerase Chain Reaction/Sequencing

[KIT Mutations in AML by Fragment Analysis and Sequencing 2002437](#)

Method: Capillary Electrophoresis

[FLT3 ITD and TKD Mutation Detection 3001161](#)

Method: Capillary Electrophoresis

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology. 500 Chipeta Way, Salt Lake City, UT 84108
(800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com
Content Review February 2021 | Last Update February 2021