

# Primary Ciliary Dyskinesia Panel

Last Literature Review: May 2021    Last Update: October 2023

Primary ciliary dyskinesia (PCD), also known as Kartagener syndrome, is a rare inherited condition that results from an underlying defect in the structure or function of motile cilia, impacting multiple body systems. Patients with PCD typically first present with neonatal respiratory distress, chronic oto-sinopulmonary disease, and year-round wet coughing. Chronic airway infection in patients with PCD often leads to bronchiectasis and obstructive lung disease, and hearing loss (either transient or permanent) can result from recurrent ear infection. Approximately half of patients with PCD will have a laterality defect such as situs inversus totalis or heterotaxy (refer to the [Heterotaxy and Situs Inversus Panel](#) Test Fact Sheet). PCD is also associated with infertility and ectopic pregnancy due to ciliary dysfunction.

## Disease Overview

### Prevalence

Approximately 1 in 16,000

### Genetics

### Etiology

Pathogenic variants in genes related to structure and function of the motile cilia

### Penetrance

100%

### Inheritance

Autosomal recessive; rare X-linked recessive forms have been reported

## Test Interpretation

### Clinical Sensitivity

Variable; dependent on phenotype/condition

### Analytic Sensitivity

For massively parallel sequencing:

| Variant Class | Analytic Sensitivity (PPA) Estimate <sup>a</sup> (%) | Analytic Sensitivity (PPA) 95% Credibility Region <sup>a</sup> (%) |
|---------------|--|--|
| SNVs          | 99.2   | 96.9-99.4  |

<sup>a</sup>Genes included on this test are a subset of a larger methods-based validation from which the PPA values are derived.

bp, base pairs; PPA, positive percent agreement; SNVs, single nucleotide variants

## Featured ARUP Testing

### [Primary Ciliary Dyskinesia Panel, Sequencing 3001621](#)

**Method:** Massively Parallel Sequencing

- Use to detect variants in genes related to primary ciliary dyskinesia, which can present with laterality defects such as situs inversus or heterotaxy as well as pulmonary disease.
- Regions of low coverage and reported variants are confirmed by Sanger sequencing as necessary.

| Variant Class       | Analytic Sensitivity (PPA) Estimate <sup>a</sup> (%) | Analytic Sensitivity (PPA) 95% Credibility Region <sup>a</sup> (%) |
|---------------------|--|--|
| Deletions 1-10 bp   | 93.8   | 84.3-98.2  |
| Deletions 11-44 bp  | 99.9   | 87.8-100   |
| Insertions 1-10 bp  | 94.8   | 86.8-98.5  |
| Insertions 11-23 bp | 99.9   | 62.1-100   |

<sup>a</sup>Genes included on this test are a subset of a larger methods-based validation from which the PPA values are derived.

bp, base pairs; PPA, positive percent agreement; SNVs, single nucleotide variants

## Limitations

- A negative result does not exclude a diagnosis of primary ciliary dyskinesia.
- Diagnostic errors can occur due to rare sequence variations.
- Interpretation of this test result may be impacted if this patient has had an allogeneic stem cell transplantation.
- The following will not be evaluated:
  - Variants outside the coding regions and intron-exon boundaries of targeted genes
  - Regulatory region and deep intronic variants
  - Noncoding transcripts
  - The following exons are not sequenced due to technical limitations of the assay:
    - ARMC4(NM\_001290020) exon 9
    - ARMC4(NM\_001290021) exon 13
    - ARMC4(NM\_001312689) exon 4
    - ARMC4(NM\_018076) exon 9
    - CCDC103(NM\_001258397) exon 4
    - CCDC114(NM\_001364171) exon 3
    - CCDC114(NM\_001364171) partial exon 4(Chr19:48822049-48822069)
    - CCDC40(NM\_001243342) exon 18
    - CFAP298(NM\_001350335) partial exon 5(Chr21:33975399-33975450)
    - CFAP298(NM\_001350337) partial exon 6(Chr21:33974534-33974561)
    - DNAAF5(NM\_017802) exon 1
    - DNAI2(NM\_001353167) exon 13
    - SPAG1(NM\_001374321) partial exon 11(Chr8:101225456-101225529)
    - SPAG1(NM\_003114) partial exon 11(Chr8:101225456-101225529)
    - SPAG1(NM\_172218) partial exon 11(Chr8:101225456-101225529)
- The following may not be detected:
  - Deletions/duplications/insertions of any size by massively parallel sequencing
  - Some variants due to technical limitations in the presence of pseudogenes, repetitive, or homologous regions
  - Low-level somatic variants

## Genes Tested

| Gene Name      | MIM #  | Associated Disorder | Inheritance Pattern |
|----------------|--------|---------------------|---------------------|
| <i>ARMC4</i>   | 615408 | PCD                 | AR                  |
| <i>CCDC103</i> | 614677 | PCD                 | AR                  |
| <i>CCDC114</i> | 615038 | PCD                 | AR                  |
| <i>CCDC151</i> | 615956 | PCD                 | AR                  |
| <i>CCDC39</i>  | 613798 | PCD                 | AR                  |

AR, autosomal recessive; XLR, X-linked recessive

| Gene Name      | MIM #  | Associated Disorder | Inheritance Pattern |
|----------------|--------|---------------------|---------------------|
| <i>CCDC40</i>  | 613799 | PCD                 | AR                  |
| <i>CCDC65</i>  | 611088 | PCD                 | AR                  |
| <i>CCNO</i>    | 607752 | PCD                 | AR                  |
| <i>CFAP298</i> | 615494 | PCD                 | AR                  |
| <i>DNAAF1</i>  | 613190 | PCD                 | AR                  |
| <i>DNAAF2</i>  | 612517 | PCD                 | AR                  |
| <i>DNAAF3</i>  | 614566 | PCD                 | AR                  |
| <i>DNAAF4</i>  | 608706 | PCD                 | AR                  |
| <i>DNAAF5</i>  | 614864 | PCD                 | AR                  |
| <i>DNAH1</i>   | 603332 | PCD                 | AR                  |
| <i>DNAH11</i>  | 603339 | PCD                 | AR                  |
| <i>DNAH5</i>   | 603335 | PCD                 | AR                  |
| <i>DNAI1</i>   | 604366 | PCD                 | AR                  |
| <i>DNAI2</i>   | 605483 | PCD                 | AR                  |
| <i>DNAL1</i>   | 610062 | PCD                 | AR                  |
| <i>DRC1</i>    | 615288 | PCD                 | AR                  |
| <i>GAS8</i>    | 605178 | PCD                 | AR                  |
| <i>LRRC6</i>   | 614930 | PCD                 | AR                  |
| <i>MCIDAS</i>  | 614086 | PCD                 | AR                  |
| <i>NME8</i>    | 607421 | PCD                 | AR                  |
| <i>PIH1D3</i>  | 300933 | PCD                 | XLR                 |
| <i>RSPH1</i>   | 609314 | PCD                 | AR                  |
| <i>RSPH3</i>   | 615876 | PCD                 | AR                  |
| <i>RSPH4A</i>  | 612647 | PCD                 | AR                  |
| <i>RSPH9</i>   | 612648 | PCD                 | AR                  |
| <i>SPAG1</i>   | 603395 | PCD                 | AR                  |
| <i>ZMYND10</i> | 607070 | PCD                 | AR                  |

AR, autosomal recessive; XLR, X-linked recessive

## Related Information

[Heterotaxy and Situs Inversus Panel](#)

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology. 500 Chipeta Way, Salt Lake City, UT 84108  
(800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com